

World In Motion


Phil Noonan
March, 2018




Overview




Realsense



Fusion




Optical Flow




OpenPose



Realsense



Fusion




Optical Flow




OpenPose

Kinect v2

**Depth Tech****Time of Flight****Depth Field of View****70.6 x 60 deg****Depth Range****0.5 m – 18.0 m
0.05 m – 0.5 m (mod)****Depth Image Resolution****512 x 424 @ 30 Hz**

R415

**Active IR Stereo****69.4 x 42.5 deg**

R435

**Active IR Stereo****91.2 x 65.5 deg****0.15 m – 10.0 m****0.2 m – 10.0 m****1280 x 720 @ 30 Hz
848 x 480 @ 90 Hz****1280 x 720 @ 30 Hz
848 x 480 @ 90 Hz**

Internals

Unit can be purchased as whole or as components

MRI compatible? Stiffener contains
stainless steel AISI 304

Magnetic Permeability 1.008 @ RT

Component	Wt. %
C	Max 0.08
Cr	18 - 20
Fe	66.345 - 74
Mn	Max 2
Ni	8 - 10.5
P	Max 0.045
S	Max 0.03
Si	Max 1


Depth Module

£58.30



Vision Processor

£25.85



R435

£155.42



R415


£128.56



Hands on, pass around!


(don't touch the lenses)

Warmup time




Static head phantom
Distance 20 cm

Warmup time



Static head phantom
Distance 25 cm
IR auto exposure on



Static head phantom
Distance 25 cm
IR auto exposure on

Both RS cameras have passive cooling


Dependent on ambient room temp

Both R415 and Kinect v2 stabilize after 20 – 30 min

R415 has varying depth of 0.6 mm vs 20+ mm for the Kinect v2 during warmup period


More investigations needed on the R435...

30 vs 90 Frames per Second




Screen captured using Nvidia shadowplay at 60 fps

30 vs 90 Frames per Second



Videos slowed down to illustrate fast twitch motion

Nyquist theorem applies!



Interference




Every 15 minutes 1 minute
of destructive interference


Light from one Kinect is
erroneously detected by the
other

Careful setup can minimize
this effect

Multipath reflections are a
big issue for ToF RGB-D






Range Modifications



Near mode – IR filter, lens refocus



Far mode – remove IR disperser,
replace wide angle lens for zoom lens






Move a slider

Other Use Scenarios

Under water

Phantoms with non-rigid internal spheres



Fish PET

Salmon fMRI

Browning et al. *Using PET/CT imaging to characterize 18 F-fluorodeoxyglucose utilization in fish*. Journal of Fish Diseases 2013

Using Mirrors

Multiple viewpoints of objects obtained


Enable views around head coil obstructions

Outdoor use

Medical imaging applications currently unclear


Plant PET?

Converse et al. *Mathematical modeling of positron emission tomography (PET) data to assess radiofluoride transport in living plants following petiolar administration*. Plant Methods 2015






Realsense



Fusion



Optical Flow



OpenPose

Iterative Closest Point







Combine sequentially obtained noisy 2.5D images into a global 3D volume

Motion between frames is small @ 30 Hz

Only few iterations needed to estimate rigid body type motion

Depth to Point

Input Data Type	Process	Output Data Type
 2D image Image pixel = float(z)	Back projection	 2D image Image pixel = float3(x, y, z)

 3D Volume Volume voxel = float(TSDF)	Raycasting	 2D image Image pixel = float3(x, y, z)
------------------------------------------------------------------------------------------------------------------------------	------------	----------------------------------------------------------------------------------------------------------------------------------

Registration

ICP is performed quickly using corresponding pixels in both float3 images


Integration

If residual error < threshold, the new depth frame is integrated into the 3D TSDF volume

Integration can be suspended to prevent small errors accumulating and causing drift



Depth to Volume



-0.9	-0.4	-0.1	0.2	0.9	1	1	1	1	1
-1	-0.5	-0.2	0.1	0.5	0.9	1	1	1	1
-1	-0.9	-0.3	0.2	0.2	0.8	1	1	1	1
-1	-0.9	-0.4	0.2	0.2	0.8	1	1	1	1
-1	-1	-0.8	-0.1	0.2	0.6	0.8	1	1	1
-1	-0.9	-0.3	-0.6	0.3	0.7	0.9	1	1	1
-1	-0.9	-0.4	-0.1	0.3	0.8	1	1	1	1
-0.9	-0.7	-0.5	0.0	0.4	0.9	1	1	1	1
-0.1	0.0	0.0	0.1	0.4	1	1	1	1	1
1	1	1	1	1	1	1	1	1	1

Registration

The TSDF is read for each 3D vertex


(potentially up to $848 \times 480 \times 90$ texture reads per second $\sim 1.1 \text{ GB s}^{-1}$)

No need to find corresponding vertices

No need to raycast

Need to perform bilinear interpolation
from the 3D texture read

OpenGL



Widely used and supported

Developed by the Khronos Group

Industry standard for multiplatform graphics

Contains methods for hardware accelerated processes

- Interpolation
- Texture/Image read+write access
- Rendering
- Plus many more ...

New API from Khronos ...



oglFusion

Aims

Opensource

Fast

Multiplatform

Upgradeable

Current Version

OpenGL 4.3 (Released August 2012)

Linux and Windows supported

CPU or GPU supported (if 4.3 compatible)

Depth to Point

Depth to Volume

Kinect v2

Realsense R415/R435

In progress

Optimisation for lower end CPU/GPU devices



Multiple camera inputs

Integrate optical flow to investigate non-rigid tracking

Long term goals


Switch to Vulkan (enabling macOS/iOS)

Multiple Camera D2V



These will get averaged out,
either they are noise only seen
by one sensor,
or the surface is too thin


0	1	2	1	0	-1
1	0	1	0	-1	-1
2	1	0	-1	-1	-2
2	1	0	-1	-2	-3
2	1	0	-1	-1	-3
3	2	1	0	-1	-2
3	2	2	1	0	-1




No problem
here, the TSDF
is consistent as
viewed from
both sensors.



Realsense



Fusion



Optical Flow



OpenPose

Optical Flow





Dense Optical Flow



Dense Optical Flow

Deepflow2




Flownet2 – neural net



Fast Dense Optical Flow

Im_0 Im_0 PyrArray Im_0 PyrArray Gradient Im_0 PyrArray Gradient Patches



Im_1 Im_1 PyrArray



Operate coarse to fine

Compute SSD from patches to Im_1

Use Hessian to converge to optimal displacement

Use previous Pyr level to initialize next level

Patch size limits resolution


Image resolution limits speed



oglFlow



openGL re-interpretation of Dense Inverse Search



Targeting non-rigid fusion



Green dots = grid
of tracking points



oglFusion requires rigid body otherwise
model blurs due to erroneous integration

With fast dense optical flow, we can acquire
knowledge of non-rigid motion before
fusion begins


oglFlow \sim 6-10 ms per frame @ 512 x 424
resolution

oglFusion \sim 1 ms per frame


90 Hz achievable in real time



Realsense



Fusion



Optical Flow



OpenPose

Suggested: Mark Ronson - Uptown Funk ft. Bruno Mars 

Real-time Multi-Person 2D Pose Estimation Using Part Affinity Fields

Zhe Cao, Tomas Simon, Shih-En Wei, Yaser Sheikh
Carnegie Mellon University

Epilepsy


Patient is monitored during hospital stay

Large amounts of video sequences

Manual annotation required

No pose quantification generally available

Patient has given approval for use of data



Detectron and DensePose



[facebookresearch / Detectron](#)

Watch 676 ⭐ Star 12,382 Fork 1,955

Code Issues 77 Pull requests 12 Projects 0 Insights

Join GitHub today
GitHub is home to over 20 million developers working together to host and review code, manage projects, and build software together.
[Sign up](#)

Dismiss


FAIR's research platform for object detection research, implementing popular algorithms like Mask R-CNN and RetinaNet.

32 commits 1 branch 0 releases 5 contributors Apache-2.0

Branch: master New pull request Find file Clone or download

ir413 and facebook-github-bot Make output dir computation dependence on datasets config fields expl... · 1 month ago

.github Add template for github issues




Example Mask R-CNN output.

DensePose: Dense Human Pose Estimation In The Wild

Rıza Alp Güler*
INRIA-CentraleSupélec
riza.guler@inria.fr

Natalia Neverova
Facebook AI Research
nneverova@fb.com

Iasonas Kokkinos
Facebook AI Research
iasonask@fb.com



DensePose-RCNN Results


DensePose COCO Dataset




Questions?




Realsense



Fusion



Optical Flow



OpenPose

Thanks





Open Heart



Laparoscopy

Analysis



Uncorrected			
Region			
Caudate	Putamen	Whole Brain	Striatum
6.80E+03	1.85E+03	118.556	1.59E+03




wRSSQ using SRTM

Kinect Based							
Region		% Decrease	Putamen	% Decrease	Whole Brain	% Decrease	Striatum
Caudate		4.77E+03	29.90	1.24E+03	32.83	9.20E+01	22.37
						7.56E+02	52.52
Image Based							
Region		% Decrease	Putamen	% Decrease	Whole Brain	% Decrease	Striatum
Caudate		4.92E+03	27.71	1.28E+03	30.79	9.36E+01	21.04
						9.20E+02	42.21

Event by Event

Pilot data using head phantom
2 internal linesource
2 minutes PET (Siemens HiRez)
Constant, handheld motion
Step changes every 5 seconds



In collaboration with
Johan Nuysts,
Matthew Bickell,
Anna Turco

